

TRANSPORT TROLLEYS - 500 SERIES

566 DN

DIMENSIONS 900 x 630 mm

HEIGHT 1310 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.
- Interior shelves for instrument boxes of 600 x 300 mm.

2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- DN Light grey.

LOAD CAPACITY

- 150 kg.

STERILISATION

- Non-airtight doors, the interior of the trolley is not isolated from contact with the outside air.

HORIZONTAL HANDLES

- 2 left side, 2 right side.

EXTENSIVE MOBILITY TROLLEYS

- Elevated push handle.
 - Dual lock wheels Ø 125 mm, 2 antistatic, 2 with brake.
- Pedal position:  Fixed wheel.  Swivel wheel.  Locked wheel.

