

CATHETER TROLLEY - 700 SERIES **720 TT.62CAT.(2)62.(4)72.73**

DIMENSIONS 900 x 630 mm

HEIGHT 1610 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

3 EXTRACTABLE SUPPORTS

- Adjustable shelf with 3 extractable supports for catheters (15 hooks).
- Identification cards ISO norm 7810, international standard.

ISO PANELS

- Flexible organisation inside the trolley.

2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- DN Light grey.
- TT Transparent.

2 VERTICAL HANDLES

- Left side.

STATIC TROLLEYS

- Double wheels Ø 100 mm, 2 antistatic, 2 with brake.



720 TT.62CAT.(2)62.(4)72.73

ACCESSORIES

62CAT CATHETERS SHELF

(2)62 ADJUSTABLE SHELVES

(4)72 ISO DRAWERS 100 mm

73 ISO DRAWER 200 mm