

STORAGE TROLLEY - 700 SERIES **720 DN.(3)62.(6)72.(2)73**

DIMENSIONS 900 x 630 mm

HEIGHT 1610 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

ISO PANELS

- Flexible organisation inside the trolley.
- Adjustable shelves, phenolic, for bulky material.

ISO DRAWERS 600 x 400 mm

- Heights 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- DN Light grey.
- TT Transparent.

2 VERTICAL HANDLES

- Left side.

STATIC TROLLEYS

- Double wheels Ø 100 mm, 2 antistatic, 2 with brake.



720 DN.(3)62.(6)72.(2)73

ACCESSORIES

(3)62 ADJUSTABLE SHELVES

(6)72 ISO DRAWERS 100 mm

(2)73 ISO DRAWERS 200 mm