

## UNIT DOSE MEDICATION - SIMPLE TROLLEY 100 SERIES 1891 B - 15 PATIENTS

**DIMENSIONS** 640 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 1005 mm (work surface – floor).

### WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

### ISO DRAWER FRAMEWORK



- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Interior frame ISO 600 x 400 mm.
- Polyethylene drawer front part, thickness 10 mm, 6 colours (Ref. G, Y, W, B, F, R).
- Ergonomic handle.
- Card holder with white cards, removable.

### ISO DRAWERS 600 x 400 mm

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

### INTERCHANGEABLE CASSETTES, 5 PATIENTS

- Compatible with pharmacy ISO furniture
- Polypropylene and aluminium, grey colour
- Cassette lock, epoxy painted steel
- Locked cassette, card upper part hidden, lower part visible.
- Cassette total weight 4,6 kg.



- Capacity: 3 litres.
- ISO norm 7810 identification card, international standard
- Appropriate for intravenous therapy medication, re-packaged medication in unit doses, pills, prefilled syringes, syrups, etc.

### KEY LOCKING

- Total locking of ISO drawers and cassettes for medication.
- Prevents from opening ISO drawers, the medication bins and the cassettes extraction.

### EXTENSIVE MOBILITY TROLLEYS

- Elevated push handle. Ergonomics, height 1045 mm.
- Double wheels Ø 125 mm, 2 antistatic, 2 with brake.



## 1891 B - 15 PATIENTS

110M6TX B.26M(R75M).(2)32L.32R.50.58

### ACCESSORIES

26M(R75M) VESA MONITOR SUPPORT

(2)32L SIDE SUPPORTS

32R SIDE SUPPORT

50 13 LITRES BIN

58 HANGER SUPPORT

### CABLE LEAD THROUGH



### KEYBOARD DRAWER

