

MRI - STORAGE TROLLEY - 700 SERIES

Trolley manufactured completely with antimagnetic materials.

706MRI OP.(2)60.(4)72.73

DIMENSIONS 640 x 480 mm

HEIGHT 1610 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

ISO PANELS

- Flexible organisation inside the trolley.
- Adjustable ISO trays, phenolic, for bulky material.

ISO DRAWERS 600 x 400 mm

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

VERTICAL HANDLES

- 2 left side.

STATIC TROLLEYS

- Double wheels Ø 100 mm, 4 antistatic, 4 with brake.



706MRI OP.(2)60.(4)72.73

ACCESSORIES

(2)60 ADJUSTABLE ISO SHELVES

(4)72 ISO DRAWERS 100 mm

73 ISO DRAWER 200 mm