

## MULTIFUNCTION - SQUARE TROLLEY SQS SERIES Q054 B

**DIMENSIONS** 480 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 980 mm (work surface – floor).

### WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

### SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

### STATIC TROLLEYS

- Without push handle - Easy left panel cleaning.
  - Free left side for mounting accessories.
- Double wheels Ø 100 mm, 2 antistatic, 2 with brake.

**Q054 B**  
**SQ061 B.32R**

### ACCESSORIES

**32R SIDE SUPPORT**

