

MULTIFUNCTION - SQUARE TROLLEY SQS SERIES Q092 W.(6)CP

DIMENSIONS 480 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 980 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

STATIC TROLLEYS

- Without push handle - Easy left panel cleaning.
 - Free left side for mounting accessories.
- Double wheels Ø 100 mm, 2 antistatic, 2 with brake.



Q092 W.(6)CP

SQ061 W.26(RA5Q.A5.R1Q).30.33L.(6)CP

ACCESSORIES

26(RA5Q.A5.R1Q) BINS DISPENSERS - CENTRAL SUPPORT

30 FOLD DOWN TRAY

33L DISPOSABLE CONTAINER SUPPORT

(6)CP PINK CARDS