

EMERGENCY - SQUARE TROLLEY SQM SERIES

Q260 R

DIMENSIONS 480 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 1015 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle. Ergonomics, height 1055 mm.
- Simple wheels Ø 125 mm, 2 antistatic, 2 with brake.



Q260 R

SQ161 R.12.14.28.32R.34.36.40

ACCESSORIES

12 MANUAL LOCKING

14 NUMBERED DISPOSABLE SEALS - 100 units

28 ROTATING ARM

32R SIDE SUPPORT

34 CARDIAC MASSAGE BOARD

36 GAS BOTTLE STAND

40 SERUM HOLDER