

TRANSPORT TROLLEYS WITH CASSETTES - 600 SERIES **608FAR.8MT.(6)C5.(2)72**

DIMENSIONS 640 x 480 mm

HEIGHT 1315 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

EPOXY PAINTED STEEL GUIDES

- Front loading, 1 side.
- C5 cassettes for transporting medications.
- 600 x 400 mm ISO drawers, heights 100, 200 mm.
- Maximum capacity 8 cassettes (medication for 40 patients).

VERTICAL HANDLES

- 2 left side, 2 right side.

EXTENSIVE MOBILITY TROLLEYS

- 4 wheel bumpers.
- Simple wheels Ø 125 mm, 2 antistatic, 2 with brake.



608FAR.8MT.(6)C5.(2)72

ACCESSORIES

8MT KEY LOCKING

(6)C5 CASSETTES 5 PATIENTS

(2)72 ISO DRAWERS 100 mm