

MRI - MULTIFUNCTION - SIMPLE TROLLEY

Trolley manufactured completely with antimagnetic materials.

305MRI

DIMENSIONS 640 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 980 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.
- Raised sides make it possible to move the trolley without needing the push handle.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

STATIC TROLLEYS

- Without push handle.
- Double wheels Ø 100 mm. Soft tread, silent.
- 4 with brake and antistatic (electrical conductive).

