

## **TREATMENT - SQUARE TROLLEY SQM SERIES**

### **Q105 B**

**DIMENSIONS** 480 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 1015 mm (work surface – floor).

#### **WORK SURFACE**

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

#### **MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Union with rivets and TIG welding.

#### **SQUARE DRAWERS**

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

#### **EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle. Ergonomics, height 1055 mm.
- Simple wheels Ø 125 mm, 2 antistatic, 2 with brake.

**Q105 B**  
**SQ132 B.32R.50**

#### **ACCESSORIES**

**32R SIDE SUPPORT**

**50 13 LITRES BIN**

