

TREATMENT - SQUARE TROLLEY SQM SERIES

Q136 B

DIMENSIONS 480 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 1015 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle. Ergonomics, height 1055 mm.
- Simple wheels Ø 125 mm, 2 antistatic, 2 with brake.



Q136 B

SQ161 B.26LW(R1Q).33R.46.50

ACCESSORIES

26LW(R1Q) CENTRAL SUPPORT

33R DISPOSABLE CONTAINER SUPPORT

46 3 GLOVES BOX HOLDER

50 13 LITRES BIN