

TREATMENT - SQUARE TROLLEY SQM SERIES

Q190 W

DIMENSIONS 480 x 480 mm

ADEQUATE HEIGHT, ERGONOMICS 1015 mm (work surface – floor).

WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Union with rivets and TIG welding.

SQUARE DRAWERS

- Polypropylene, PP. 3 colours - B Blue, W White, R Red.
- Heights 57, 130 mm.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 20 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle. Ergonomics, height 1055 mm.
- Simple wheels Ø 125 mm, 2 antistatic, 2 with brake.



Q190 W

SQ161TXQ W.26M(R75QM).30.32L.46.50

ACCESSORIES

26M(R75QM) VESA MONITOR SUPPORT TILT AN ROTATION ADJUSTABLE

30 FOLD DOWN TRAY

32L SIDE SUPPORT

46 3 GLOVES BOX HOLDER

50 13 LITRES BIN

CABLE LEAD THROUGH



KEYBOARD DRAWER

