

TRANSPORT TROLLEYS - TROLLEY WITH FIXED SHELVES 603 DN

DIMENSIONS 640 x 480 mm
HEIGHT 1655 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.
- 2 central fixed shelves.

2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys.
- Card holder with white cards, removable, left door.
- **DN** Light grey.

LOAD CAPACITY

- 150 kg.

VERTICAL HANDLES

- 2 left side, 2 right side.

EXTENSIVE MOBILITY TROLLEYS

- 4 wheel bumpers for hits.
 - Wheels load capacity 400 kg.
 - Ø 125 mm double wheels, 2 antistatic, 2 with dual lock.
- Pedal position:  Fixed wheel.  Swivel wheel.  Locked wheel.

