

## TRANSPORT TROLLEYS - TROLLEY WITH FIXED SHELVES 603 OP

**DIMENSIONS** 640 x 480 mm  
**HEIGHT** 1655 mm

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.
- 2 central fixed shelves.

### LOAD CAPACITY

- 150 kg.

### VERTICAL HANDLES

- 2 left side, 2 right side.

### EXTENSIVE MOBILITY TROLLEYS

- 4 wheel bumpers for hits.
  - Wheels load capacity 400 kg.
  - Ø 125 mm double wheels, 2 antistatic, 2 with dual lock.
- Pedal position:  Fixed wheel.  Swivel wheel.  Locked wheel.

