

TRANSPORT TROLLEYS - OPEN TROLLEY **610TRT.(8)72.73**

DIMENSIONS 640 x 480 mm

HEIGHT 1650 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Assembly with rivets and TIG welding.

EPOXY PAINTED STEEL GUIDES

- Front loading, 1 side.

ISO DRAWERS 600 x 400 mm

- Heights 50, 100, 200 mm. Interchangeable, ABS.
- ISO norm 3394.

LOAD CAPACITY


- 100 kg.

VERTICAL HANDLES

- 2 left side, 2 right side.

EXTENSIVE MOBILITY TROLLEYS

- 4 wheel bumpers for hits.
- Ø 125 mm simple wheels, 2 antistatic, 2 with brake.



610TRT.(8)72.73

ACCESSORIES

(8)72 ISO DRAWERS 100 mm

73 ISO DRAWER 200 mm