

1

2

3

4

A


SERUM SIDE SUPPORT / SOPORTE LATERAL PARA SUERO

SQUARE TROLLEY / CARRO CUADRADO 480 x 480 mm**SIMPLE TROLLEY / CARRO SIMPLE 640 x 480 mm**

- REF. 32SL - LEFT SIDE / LADO IZQUIERDO
- REF. 32SR - RIGHT SIDE / LADO DERECHO

- MATERIAL: GREY EPOXY PAINTED STEEL
- 2 FIXATION SLEEVES IN THE TROLLEY
- MANUALLY REMOVABLE
- EASY-TO-CLEAN TROLLEY STAINLESS STEEL PANEL
- LOAD CAPACITY: 8 SERUM BOTTLES OF 1L.

- MATERIAL: ACERO PINTADO EPOXI GRIS
- 2 MANGUITOS DE FIJACIÓN EN EL CARRO
- DESMONTABLE MANUALMENTE
- FÁCIL LIMPIEZA DEL PANEL DE ACERO INOXIDABLE DEL CARRO
- CAPACIDAD DE CARGA: 8 BOTELLAS DE SUERO DE 1L.



B

C

D



06		Dibujado	Fecha	ESCALA	 INSAUSTI®	
05			24-09-25	1:4		
04		Comprobado	Fecha			
03						
02						
01	14-11-25	FICHA TÉCNICA REF. 32SL				PLANO Nº
MODIFICACIONES						350-0324-01